


# Quick Reference Guide Masimo Root VSM

DISCLAIMER: This document has not been approved by Masimo Marketing and was created for your site use ONLY. Not to be shared with other healthcare facilities

1. Place SpO2 sensor on non-dominant ring or middle finger



2. Click on the top left corner of the screen to start monitoring



OR

Push the power button on the front of the handheld monitor and push the side button of the Root monitor



3. Changing Alarm Parameters



- Click on the parameter you want to change (ex: SpO2)
- Click "Alarms"
- Change high and low limit → select OK after change

4. Blood Pressure



- Click the icon of the patient with cuff illustration
- When illuminated in blue, the cuff is inflating

## To set NIBP Intervals



- Tap drop down button at top right corner of NIBP box
- Tap Intervals
- Select desired interval time
- Press OK
- Press Start Auto

**NOTE:** BP cuff must be left on patient's arm for interval measurement

5. Display

- To change views click drop down arrow indicated below
- Click trend view to review past vital signs data



6. Remove SpO2 sensor

- Select STANDBY button top left-hand corner, select confirm

