

SMARTGUARD™ AUTO MODE TRAINING CHECKLIST MINIMED™ 770G SYSTEM

Medtronic

Certified Product Trainer: _____

Initial Pump Training Date: _____ AM Training Date: _____ Overall time spent: _____

Patient Name: _____

Pump S/N#: _____ Completed Online Training Video(s) Yes No

REMINDERS FOR SUCCESS HAVE BEEN REVIEWED:

- Regular fingerstick BGs and calibrations are needed; if BG result seems inaccurate, wash hands and re-check
- Give boluses before all meals and snacks; give correction boluses when the pump recommends
- Perform bedtime calibration and give correction bolus if needed to minimize nighttime alerts
- If SmartGuard™ Auto Mode exit occurs, the goal is to re-enter as soon as possible
- Regular CareLink™ Personal software uploads are needed

SmartGuard™ Auto Mode

- Introduction to SmartGuard™ Auto Mode compared to Manual Mode

CareLink™ Software – Manual Mode Evaluation

- Calibrations are entered appropriately
- Carbs are being accurately quantified
- Boluses given pre-meal/snack
- Post-meal BGs within target
- Correction boluses achieving target
- Basal/Bolus ratio evaluated

Actions taken if any of these not achieved: _____

- CareLink™ Software report is not available

Ensure SmartGuard™ Auto Mode Readiness

- SmartGuard™ Auto Mode warm-up is complete
- Carb Ratio and Active Insulin time entered
- A sensor is on and working
- There is no active Temp Basal, Manual Suspend or Square/Dual Wave bolus

SmartGuard™ Auto Mode Start

- SmartGuard™ Auto Mode turned On
- SmartGuard™ Auto Mode Readiness screen checked
- BG has been entered
- SmartGuard™ Auto Mode BG alert is On
- All settings entered as per the MiniMed™ 770G system initiation settings form

Additional topics

- SmartGuard™ Auto Mode Home screen/graphs
- Safe Basal
- SmartGuard™ Auto Mode exit reasons reviewed
- Calibration
- Follow instructions on screen to resolve
- Check SmartGuard™ Auto Mode Readiness screen for re-entry information
- Reviewed common alerts
- Temp Target
- Bolusing in SmartGuard™ Auto Mode

Resources

- Getting Started with MiniMed™ 770G SmartGuard™ Auto Mode
- Follow-up touch points scheduled

SmartGuard™ Auto Mode has been activated: Yes No If no, SmartGuard™ Auto Mode start scheduled for: _____

Training Notes: _____

I certify that I received and understood the full and complete training mentioned above

Patient Signature: _____

Patient Name: _____

Follow-up Touch Point Date: _____

Assessment topic and observations:

% sensor wear _____ % time in SmartGuard™ Auto mode _____ Carb Ratio: Time Ratio

% time in range _____ % auto basal _____ _____ _____

SmartGuard™ Auto Mode exits Bolus timing Correction boluses _____ _____

Temp Target Overnight alerts Active insulin time _____ _____ _____

Manual Mode setting review _____ _____

Action taken to address the observations above: _____

Follow-up Touch Point Date: _____ Active Insulin Time _____

Assessment topic and observations:

% time in SmartGuard™ Auto Mode _____ % Sensor Wear: _____ Carb Ratio: Time Ratio

Glucose control when in Manual Mode Meal boluses Correction boluses _____ _____

Night-time alerts Activity/exercise SmartGuard™ Auto Mode exits _____ _____

Action taken to address the observations above: _____

Follow-up Touch Point Date: _____ Active Insulin Time _____

Assessment topic and observations:

% time in SmartGuard™ Auto Mode _____ % Sensor Wear: _____ Carb Ratio: Time Ratio

Glucose control when in Manual Mode Meal boluses Correction boluses _____ _____

Night-time alerts Activity/exercise SmartGuard™ Auto Mode exits _____ _____

Action taken to address the observations above: _____

By signing below the patient acknowledges that they have received and understood the full and complete training mentioned above and consents to a copy of this document being forwarded to Medtronic for administrative and quality purposes only.

Patient Signature: _____ Date: _____

Trainer signature: _____ Date: _____

FAX COMPLETED FORM WITHIN 30 DAYS OF TRAINING TO 1-800-326-8322 or send scan/photo (in pdf format) to your specific Region's email address: rs.cancptinvoiceon@medtronic.com for ON, rs.can-cptinvoice-west@medtronic.com for BC/AB/MB/SK, rs.can-cptinvoice-east@medtronic.com for NS/NB/NL/PEI, rs.can-cptinvoice-qc@medtronic.com for QC