

MiniMed™ 780G system detailed menu map



History & Graph

History
Sensor Glucose Review
Graph
Time in Range



SmartGuard

SmartGuard Checklist
Temp Target
SmartGuard Settings
SmartGuard On and Off



Sound & Vibration

Silence Sensor Alerts
Volume
Sound
Vibration
Alert Settings Shortcut



Reservoir & Set

New Reservoir
Fill Cannula



Insulin

Bolus
Basal
Suspend/Resume Basal
Delivery
Delivery Settings Shortcut



Blood Glucose

BG



Status

Suspend All Delivery
SmartGuard Checklist
Pump
Sensor



Paired devices

Pair New Device
Pair CareLink
Sensor



Settings

Alert Settings

- High Alert
- Low Alert
- Snooze High & Low
- Reminders

Delivery Settings

- Bolus Wizard Setup
- Basal Pattern Setup
- Max Basal/Bolus
- Dual/Square Wave
- Bolus Increment
- Bolus Speed
- Preset Bolus Setup
- Preset Temp Setup

Device Settings

- Sensor On and Off
- Time & Date
- Device Info
- Display
- Block Mode
- Self Test
- Review Settings
- Manage Settings
- Easy Bolus
- Auto Suspend
- Language